

# A-dec Choices

Which solution is right for you?



a dec<sup>®</sup>



## A-dec 500 | Premium. Graceful. Comfort.

**When you need the best.** A-dec 500 is dentistry's hallmark for treatment room comfort *and* performance. A-dec 500 gives you optimal ergonomics supported by a full range of integrated possibilities. And with a chair ride that is exquisitely smooth, you will see why A-dec 500 is unmatched in its graceful innovation.

- Outstanding comfort and support, plus an ultra-thin 25mm backrest for optimal oral cavity access and increased leg room.
- Seamless ancillary integration to meet the changing needs of dentistry.
- Anatomically shaped contours reduce patient pressure points.
- Exceptional range of motion plus low "down" position improves access, visibility and posture.
- Integrated touchpad to give you easy control of chair, cuspidor, and dental light.
- Control head pivot points provide superb instrument positioning.

## A-dec 400 | Accessible. Strong. Style.

**Contemporary design in a versatile package.** Modest luxury with proven reliability. A-dec 400 marries form and function for superb patient and doctor comfort and access. The ambidextrous compatibility of the Radius® delivery system brings flexibility and productivity to every treatment room. And when you consider its unique load-bearing pivoting armrests, contemporary chair design, and 60° swivel, you will see why A-dec 400 is an exceptional value.

- Ultra-thin backrest enables you to tuck your knees under the chair for optimal ergonomic access.
- Robust two-position armrest design gives secure patient support, and stows easily out of the way for direct access by the dental team.
- Vertical range of 349-800mm accommodates stand-up dentistry.
- Radius delivery and support modules rotate quickly and easily around the chair for complete left/right compatibility.

## A-dec 300 | Smart. Sensible. Flexible.

**Proven performance for any budget.** Start with the chair: optional armrests and sleek seamless or plush sewn upholstery. Choose a Traditional or Continental delivery system to suit your working style. With Radius left/right assistant's arm, a wide range of ancillary possibilities, and options for touchpads and foot controls, you can truly create a smart system that fits your practice—and your budget.

- Simplified design with fewer moving parts reduces maintenance and extends product life.
- Ultra-thin 32mm backrest means good oral cavity access and improved posture.
- Radius delivery option for versatile left/right compatibility.
- Smooth, quiet hydraulic lift cylinder engineered for a 20 year life.
- Seamless ancillary integration allows additional equipment as your needs change
- Modular system can grow and evolve right along with your practice
- Vertical range of 349-800mm accommodates stand-up dentistry.

## A-dec 200 | Proven. Efficient. Smooth.

**Reliability in a complete package.** Genuine A-dec for any budget. Designed and engineered to withstand years of constant use, A-dec 200 is a technology-ready solution that includes dental chair, delivery system, assistant's instrumentation, dental light, and support center with cuspidor.

- Double-articulating headrest with easy adjustment.
- Locking 2-position swiveling armrests for easy patient entry and exit.
- Reliable hydraulic-lift system for ultra-quiet, smooth chair movement.
- 300-watt power supply.
- Optional total touchpad control.

<b>CHAIR</b>	<b>A-dec 200</b>	<b>A-dec 300</b>	<b>A-dec 400</b>	<b>A-dec 500</b>
Choice of sewn or seamless upholstery	Seamless only	•	•	•
Dual-articulating gliding headrest	Locking Knob (without glide)	Locking Knob	Locking Knob or Lever Release	Lever Release
Thin backrest (25-32mm)		•	•	•
Virtual pivot for added patient comfort		•	•	•
Ultra-soft start and stop drive system				•
Multi-position armrest settings	2	2	2	3
Swivel (30° each side of center)			•	•
Hydraulic drive	•			•
Hybrid drive (hydraulic base/electric tilt)		•	•	
Chair at lowest point	457mm	349mm	349mm	343mm
Chair at highest point	737mm	800mm	800mm	800mm
Integrated 300W system power supply	•	•	•	•
<b>DELIVERY SYSTEM</b>				
Dental system modules	A-dec 200 A-dec 300	A-dec 300	A-dec 300 A-dec 400	A-dec 300 A-dec 400 A-dec 500
Three-position control block	Optional	Standard	Standard	—
Four-position control block	Standard	Optional	Optional	Standard
Ancillary accommodation	Yes	2	2	4
Choice of programmable standard or deluxe touchpad with DCS	•	•	•	•
Self-contained water system with QD 2-liter bottle	Threaded 2-liter (no QD)	•	•	•
Accommodates amalgam separator	•	•	•	•
Cuspidor swivel	•		•	•
Programmable timing for cup and bowl	Optional	•	•	•
<b>DENTAL LIGHT</b>	<b>A-dec 200 Halogen 2-Axis</b>	<b>A-dec 500 Halogen 3-Axis</b>	<b>A-dec 300 LED LED 3-Axis</b>	<b>A-dec 500 LED LED 3-Axis</b>
Auto on/off (using chair pre-position)	•	•	•	•
Composite mode	•	•	•	Cure-safe
Intensity settings	Composite: 5,500 lux High: 20,000 lux	Composite: 8,000 lux Mid: 20,000 lux High: 24,000 lux	Composite: 8,000 lux Low: 8,000 lux High: 25,000 lux	Low: 15,000 lux Mid: 25,000 lux High: 30,000 lux Cure-safe: 25,000 yellow light
Color temperature	4,800K	5,000K	5,000K	5,000K
Power consumption	95W	95W	15W	15W



## A-dec Worldwide

With A-dec, you gain a partner for the life of your practice. From our headquarters in Newberg, Oregon, we support global offices and a network of dedicated dealers in more than 100 countries worldwide.

Discover the A-dec Difference. Visit [a-dec.com/TheDifference](http://a-dec.com/TheDifference) to find out how every detail behind an A-dec solution ensures a reliable, lasting investment.



### A-dec Headquarters

2601 Crestview Drive  
Newberg, Oregon 97132 USA  
Tel: 1.800.547.1883 within USA/CAN  
Tel: +1.503.538.7478 outside USA/CAN  
Fax: +1.503.538.0276  
[a-dec.com](http://a-dec.com)

### A-dec Australia

Tel: 1.800.225.010 within AUS  
Tel: +61.(0)2.8332.4000 outside AUS

### A-dec China

Tel: 400.600.5434 within CN  
Tel: +86.571.89026088 outside CN

### A-dec United Kingdom

Tel: 0800.233285 within UK  
Tel: +44.(0).24.7635.0901 outside UK

A-dec, the A-dec logo, Cascade, Century Plus, Continental, Decade, ICX, ICV, Performer, Preference, Preference Collection, Preference ICC, and Radius are trademarks of A-dec Inc. and are registered in the United States and other countries. A-dec 500, A-dec 400, A-dec 300, A-dec 200, A-dec Inspire, and reliablecreativesolutions are also trademarks of A-dec Inc. None of the trademarks or trade names in this document may be reproduced, copied, or manipulated in any manner without the express, written approval of the trademark owner.

Certain touchpad symbols are proprietary to A-dec Inc. Any use of these symbols, in whole or in part, without the express written consent of A-dec Inc., is strictly prohibited.

The colors and products shown may vary from the actual materials and are subject to change without notice. Please see your authorized A-dec dealer for accurate color samples and the most current product information.